

Product Datasheet

Canine Interleukin 2 ELISA Kit (orb1214681)

Description	Canine Interleukin 2, IL-2 ELISA Kit
Reactivity	Canine
Range	6-1800ng/L
Tested Applications	ELISA
Target	IL2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	17,668
Uniprot ID	Q29416
Solubility (25°C)	1920ng/L
Sensitivity	3.43ng/L
Entrez	403989
Expiration Date	Please enquire.