

Product Datasheet

Canine Endothelin 1 ELISA Kit (orb1214697)

Description	Canine Endothelin 1, ET-1 ELISA Kit
Reactivity	Canine
Range	1-400ng/L
Tested Applications	ELISA
Target	ET-1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	22,977
Uniprot ID	P13206
Solubility (25°C)	480ng/L
Sensitivity	0.53ng/L
Entrez	403424
Expiration Date	Please enquire.