

Product Datasheet

Canine Caveolin-1 ELISA Kit (orb1214699)

Description	Canine Caveolin-1, CAV-1 ELISA Kit
Reactivity	Canine
Range	20-6000pg/ml
Tested Applications	ELISA
Target	CAV1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	20,606
Uniprot ID	P33724
Solubility (25°C)	6400pg/ml
Sensitivity	10.68pg/ml
Entrez	403980
Expiration Date	Please enquire.