

Product Datasheet

Canine Interferon Gamma ELISA Kit (orb1214700)

Description	Canine Interferon Gamma, IFN-G ELISA Kit
Reactivity	Canine
Range	5-1500ng/L
Tested Applications	ELISA
Target	IFNG
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	19,395
Uniprot ID	P42161
Solubility (25°C)	1920ng/L
Sensitivity	2.58ng/L
Entrez	403801
Expiration Date	Please enquire.