

Product Datasheet

Canine Interleukin 10 ELISA Kit (orb1214701)

Description	Canine Interleukin 10, IL-10 ELISA Kit
Reactivity	Canine
Range	2-700pg/ml
Tested Applications	ELISA
Target	IL10
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	20,906
Uniprot ID	P48411
Solubility (25°C)	800pg/ml
Sensitivity	1.02pg/ml
Entrez	403628
Expiration Date	Please enquire.