

Product Datasheet

Canine Interleukin 6 ELISA Kit (orb1214702)

Description	Canine Interleukin 6, IL-6 ELISA Kit
Reactivity	Canine
Range	0.05-15ng/L
Tested Applications	ELISA
Target	IL6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	22,945
Uniprot ID	P41323
Solubility (25°C)	16ng/L
Sensitivity	0.02ng/L
Entrez	403985
Expiration Date	Please enquire.