

Product Datasheet

Bovine Tnf Receptor-associated Factor 6 ELISA Kit (orb1214717)

Description	Bovine Tnf Receptor-associated Factor 6, TRAF6 ELISA Kit
Reactivity	Bovine
Range	15-3000ng/L
Tested Applications	ELISA
Target	TRAF6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	61,647
Uniprot ID	Q3ZCC3
Solubility (25°C)	3200ng/L
Sensitivity	6.48ng/L
Entrez	539124
Expiration Date	Please enquire.