

Product Datasheet

Bovine Heat Shock Protein Beta-1 ELISA Kit (orb1214732)

Description	Bovine Heat Shock Protein Beta-1, HSPB1 ELISA Kit
Reactivity	Bovine
Range	0.05-30ng/mL
Tested Applications	ELISA
Target	HSPB1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	22,393
Uniprot ID	Q3T149
Solubility (25°C)	32ng/mL
Sensitivity	0.021ng/mL
Entrez	516099
Expiration Date	Please enquire.