

## Product Datasheet

# Bovine C-C Motif Chemokine 5 ELISA Kit (orb1214739)

<b>Description</b>	Bovine C-C Motif Chemokine 5, CCL5 ELISA Kit
<b>Reactivity</b>	Bovine
<b>Range</b>	2-600ng/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	CCL5
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	10,059
<b>Uniprot ID</b>	<a href="#">O97919</a>
<b>Solubility (25°C)</b>	640ng/mL
<b>Sensitivity</b>	1.27ng/mL
<b>Entrez</b>	<a href="#">327712</a>
<b>Expiration Date</b>	Please enquire.