

Product Datasheet

Bovine Gamma-Interferon-inducible Lysosomal Thiol Reductase ELISA Kit (orb1214768)

Description	Bovine Gamma-Interferon-inducible Lysosomal Thiol Reductase, IFI30 ELISA Kit
Reactivity	Bovine
Range	0.05-15ng/ml
Tested Applications	ELISA
Target	IFI30
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	27,490
Uniprot ID	A6QPN6
Solubility (25°C)	16ng/ml
Sensitivity	0.01ng/ml
Entrez	615930
Expiration Date	Please enquire.