

Product Datasheet

Bovine TGF-beta Receptor Type-1 ELISA Kit (orb1214783)

Description	Bovine TGF-beta Receptor Type-1, TGFBR1 ELISA Kit
Reactivity	Bovine
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	TGFBR1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	55,814
Uniprot ID	O46680
Solubility (25°C)	64ng/ml
Sensitivity	0.092ng/ml
Entrez	282382
Expiration Date	Please enquire.