

Product Datasheet

Bovine Interleukin 1B ELISA Kit (orb1214785)

Description	Bovine Interleukin 1B, IL-1B ELISA Kit
Reactivity	Bovine
Range	0.5-200ng/L
Tested Applications	ELISA
Target	IL1B
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	30,760
Uniprot ID	P09428
Solubility (25°C)	240ng/L
Sensitivity	0.25ng/L
Entrez	281251
Expiration Date	Please enquire.