

Product Datasheet

Bovine Apelin ELISA Kit (orb1214790)

Description	Bovine Apelin, AP ELISA Kit
Reactivity	Bovine
Range	5-2000ng/ml
Tested Applications	ELISA
Target	APLN
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	8,743
Uniprot ID	Q9TUI9
Solubility (25°C)	2400ng/ml
Sensitivity	2.47ng/ml
Entrez	282143
Expiration Date	Please enquire.