

Product Datasheet

Bovine Vascular Endothelial Cell Growth Factor A ELISA Kit (orb1214799)

Description	Bovine Vascular Endothelial Cell Growth Factor A, VEGFA ELISA Kit
Reactivity	Bovine
Range	20-6000ng/L
Tested Applications	ELISA
Target	VEGFA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	22,310
Uniprot ID	P15691
Solubility (25°C)	6400ng/L
Sensitivity	9.08ng/L
Entrez	281572
Expiration Date	Please enquire.