

Product Datasheet

CD80 antibody (PE/Cy7) (orb1214863)

Description	Hamster monoclonal antibody against CD80 conjugated to PE/Cy7
Species/Host	Hamster
Reactivity	Mouse
Conjugation	PE/Cy7
Tested Applications	FC
Storage	2-8°C protected from light
Note	For research use only
Isotype	Armenian Hamster IgG
Clonality	Monoclonal
Uniprot ID	Q00609
Entrez	12519
Expiration Date	12 months from date of receipt.