

## Product Datasheet

### F4/80 antibody (PE/Cy7) (orb1215160)

<b>Description</b>	Rat monoclonal antibody against F4/80 conjugated to PE/Cy7
<b>Species/Host</b>	Rat
<b>Reactivity</b>	Mouse
<b>Conjugation</b>	PE/Cy7
<b>Tested Applications</b>	FC
<b>Storage</b>	2-8°C protected from light
<b>Note</b>	For research use only
<b>Isotype</b>	Rat IgG2b, kappa
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">Q61549</a>
<b>Entrez</b>	<a href="#">13733</a>
<b>Expiration Date</b>	12 months from date of receipt.