

## Product Datasheet

### Rat IgG2a Isotype Control (PE/Cy5) (orb1215214)

<b>Description</b>	Rat monoclonal Rat IgG2a Isotype Control conjugated to PE/Cy5
<b>Species/Host</b>	Rat
<b>Conjugation</b>	PE/Cy5
<b>Tested Applications</b>	FC
<b>Storage</b>	2-8°C protected from light
<b>Note</b>	For research use only
<b>Isotype</b>	Rat IgG2a
<b>Clonality</b>	Monoclonal
<b>Expiration Date</b>	12 months from date of receipt.