

Product Datasheet

CD11b antibody (PE/Cy5) (orb1215221)

| | |
|----------------------------|--|
| Description | Rat monoclonal antibody against CD11b conjugated to PE/Cy5 |
| Species/Host | Rat |
| Reactivity | Human, Mouse |
| Conjugation | PE/Cy5 |
| Tested Applications | FC |
| Storage | 2-8°C protected from light |
| Note | For research use only |
| Isotype | Rat IgG2b, kappa |
| Clonality | Monoclonal |
| Entrez | 16409 |
| Expiration Date | 12 months from date of receipt. |