

Product Datasheet

In Vivo Anti-Mouse CD40 (FGK45) (orb1215339)

Description	In Vivo Anti-Mouse CD40 (FGK45)
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FA, FC
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rat IgG2a
Clonality	Monoclonal
Clone Number	FGK45
Antibody Type	Primary Antibody
Entrez	21939
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)