

## Product Datasheet

# Mouse IgG2a Isotype Control antibody (orb1215351)

<b>Description</b>	Mouse monoclonal to Mouse IgG2a Isotype Control
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FA, FC, IHC-Fr, IHC-P, IP
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Isotype</b>	Mouse IgG2a
<b>Clonality</b>	Monoclonal
<b>Expiration Date</b>	12 months from date of receipt.