

## Product Datasheet

### Mouse IgG2a Isotype Control (FITC) (orb1215399)

<b>Description</b>	Mouse monoclonal Mouse IgG2a Isotype Control conjugated to FITC
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	FITC
<b>Tested Applications</b>	FC
<b>Storage</b>	2-8°C protected from light
<b>Note</b>	For research use only
<b>Isotype</b>	Mouse IgG2a
<b>Clonality</b>	Monoclonal
<b>Expiration Date</b>	12 months from date of receipt.