


Product Datasheet

MPXV A35R Protein (orb1215649)

Description	Recombinant MPXV A35R Protein is produced by mammalian expression system and the target gene encodin...
Reactivity	Virus
Tested Applications	ELISA, MS, SDS-PAGE, WB
Preservatives	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.
Form/Appearance	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Please minimize freeze-thaw cycles.
Purity	Greater than 95% as determined by reducing SDS-PAGE.
Source	Human cells
Uniprot ID	AUW64225.1
Expiration Date	6 months from date of receipt.



Greater than 95% as determined by reduci...