


## Product Datasheet

### MPXV A29L Protein (orb1215650)

<b>Description</b>	Recombinant MPXV A29L Protein is produced by E.coli expression system and the target gene encoding S...
<b>Reactivity</b>	Virus
<b>Tested Applications</b>	ELISA, MS, SDS-PAGE, WB
<b>Preservatives</b>	Liquid in a 0.2 µm filtered solution of PBS, 300 mM NaCl, pH 7.4.
<b>Form/Appearance</b>	Liquid in a 0.2 µm filtered solution of PBS, 300 mM NaCl, pH 7.4.
<b>Storage</b>	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Please minimize freeze-thaw cycles.
<b>Purity</b>	Greater than 95% as determined by reducing SDS-PAGE.
<b>Source</b>	E.coli
<b>Uniprot ID</b>	<b>P48865</b>
<b>Expiration Date</b>	6 months from date of receipt.



Greater than 95% as determined by reduci...