

## Product Datasheet

### Mammal FLT3 Ligand protein (orb1215653)

<b>Description</b>	The Jamaican fruit-eating Bat FLT3 Ligand yeast-derived recombinant protein is not tagged and is nat...
<b>Target</b>	FLT3 Ligand
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	RTPDCTFPHS PISSTFADTI RKLSDYLLLED YPVTVATNLH DDKLCWAFWR LVLAQRCMER LKTVAGSQMQ QLLEDVNTEI YFVTSCALQD TSEQLVALKP WITRWNFSGC LELRCQPDSS TLLPSRSPKA LEATALPAAQ AP (142)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	16.0 kDa
<b>Entrez</b>	<b>119046007</b>
<b>Expiration Date</b>	6 months from date of receipt.