

## Product Datasheet

### Rabbit IL-25 protein (orb1215665)

<b>Description</b>	The Rabbit IL-25 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
<b>Target</b>	IL-25
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	LELWKDCSHR STCCTSKGQA LPEEWLQWTP VPVTPPELEL THPQESCRAS EDGPLNSRAI SPWRYELDRD LNRLPQDLYH ARCLCPHCVS LQTGSHMDPL GNSELLYHNQ TVFYRRPCRG ERGLRQGYCW KEDCIEFLWL VFVLDQELWL K (151)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	17.4 kDa
<b>Entrez</b>	<b>100345869</b>
<b>Expiration Date</b>	6 months from date of receipt.