

## Product Datasheet

### Bovine TIM-3 protein (orb1215674)

<b>Description</b>	The Bovine TIM-3 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
<b>Target</b>	TIM-3
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	AYVSQVGQNA DLPCTYSPAT TENLVPVCWG KGPCPVFECY SLVLRDGRN VTYQTSSRYL LKRDLHKGDV TLTIKNVTLA DSGTYCCRIQ FPGLMNRKS NLELIKPAK VTPAWTPWRD ITTAFFRMLT TKGPVSETRT KTLHDKNQT EISTLATELQ DMGATTRTG (169)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	18.8 kDa
<b>Entrez</b>	<b>768047</b>
<b>Expiration Date</b>	6 months from date of receipt.