

## Product Datasheet

### Guinea pig TNF alpha protein (orb1215677)

<b>Description</b>	The Guinea Pig TNF alpha Biotinylated yeast-derived recombinant protein is not tagged and is natural...
<b>Target</b>	TNF alpha
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	SASQNDNDKP VAHVVANQQA EEELQWLSKR ANALLANGMG LSDNQLVVP DGLYLIYSQV LFKGQGCPY LLLTHTVSRL AVSYPEKVNLSAIKSPCQK ETPEGAERKP WYEPIYLGGV FQLQKGDRLS AEVNLPQYLD FADSGQIYFG VIAL (154)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	17.0 kDa
<b>Entrez</b>	<b>100135630</b>
<b>Expiration Date</b>	6 months from date of receipt.