

Product Datasheet

Feline IL-4 protein (orb1215680)

Description	The Feline IL-4 Biotinylated yeast-derived recombinant protein is not tagged and is naturally endoto...
Target	IL-4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QNFNNTLKEI IKTLNILTAR NDSCMELTVM DVLAAPKNTS DKEIFCRATT VLRQIYTHHN CSTKFLKGLD RNLSSMANRT CSVNEVKKCT LKDFLERLKA IMQKKYSKH (109)
Purity	98%
Source	Yeast
MW	12.6 kDa
Entrez	751514
Expiration Date	6 months from date of receipt.