

## Product Datasheet

### Rabbit IL-4 protein (orb1215688)

<b>Description</b>	The Rabbit IL-4 Biotinylated yeast-derived recombinant protein is not tagged and is naturally endoto...
<b>Target</b>	IL-4
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	GRRGDILPE VIKTLNILTE RKTPTKLM I ADALAVPKNT TEREAVCRAA TALRQFYLHH KVSWC FKEHG ELGDLRLLRG LDRNLCSMAK LSNCPGKEAR QTTLEDFLDR LKTAMQEKYS KRQS (124)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	14.2 kDa
<b>Entrez</b>	<b>100302454</b>
<b>Expiration Date</b>	6 months from date of receipt.