

Product Datasheet

Porcine IL-22 protein (orb1215689)

Description	The Swine IL-22 Biotinylated yeast-derived recombinant protein is not tagged and is naturally endoto...
Target	IL-22
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	VPITHHCKLD QSNFQQPYIT NRTFTLAQEA SLADNNTDVR LIGNNLFQGV NMRERCYLVK QVLNFTLEEV LFPNSDRFHP YMQEVASFLD SLSKKLSQCR IKGDDQHIQR NVNFKDTVK KLGESGEIKV IGELYLLFMA LKNECTLPGH SWKMDN (156)
Purity	98%
Source	Yeast
MW	18.0 kDa
Entrez	595104
Expiration Date	6 months from date of receipt.