

## Product Datasheet

### Gallus IL-12 p40 protein (orb1215695)

<b>Description</b>	The Chicken IL-12/IL-23 p40 Biotinylated yeast-derived recombinant protein is not tagged and is natu...
<b>Target</b>	IL-12 p40
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	AQWKLRENVY VIESEWNDET PAKKVKLTCD TSDEALPVYW KKGTELKGTG KTLTTEVKEF PDAGNYTCLS AKTHEIISYS FFLITKVDSN GQMIRSILKS YKEPSKFLK CEAKNYSGIF TCSWMTENES PSVKFTIRSL KGSQGDVTCS SPVARTDKSV TEYTAQCQKE NYCPFAEEHQ PTEMFLEVID EVEYENYTSS FFIRDIKPD PPQCQYASTN GTVTWTYPKT WSTPKSYFPL TFRVKVESTK KYKSKVYDAD EQSIQIPKTG PKDKISVQAR DRYYNSSWSE WSTLCR (296)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	34.1 kDa
<b>Entrez</b>	<a href="#">404671</a>
<b>Expiration Date</b>	6 months from date of receipt.