

Product Datasheet

Canine TNF alpha protein (orb1215698)

Description	The Canine TNF alpha Biotinylated yeast-derived recombinant protein is not tagged and is naturally e...
Target	TNF alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	VKSSSRTPSD KPVAHVVANP EAEGQLQWLS RRANALLANG VELTDNQLIV PSDGLYLIYS QVLFKGGQCP STHVLLHTI SRFAVSYQTK VNLLSAIKSP CQRETPEGTE AKPWYEPIYL GGVFQLEKGD RLSAEINLPN YLDFAESGQV YFGIIAL (157)
Purity	98%
Source	Yeast
MW	17.3 kDa
Entrez	403922
Expiration Date	6 months from date of receipt.