

## Product Datasheet

### Rabbit CTLA-4 protein (orb1215701)

<b>Description</b>	The Rabbit sCTLA-4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. ...
<b>Target</b>	CTLA-4
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	LHVSQPAVVLASSRGVASFVCEYASSHKATEVRVTVLRQANSQMTEVCAMTYTV ENELTFIDDSTCTGISHGNKVNLTIQGLSAMDTGLYICKVELMYPPPPYYVGMGNG TQIYVIEPEPCPDS (124)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	13.5 kDa
<b>Entrez</b>	<a href="#">100009412</a>
<b>Expiration Date</b>	6 months from date of receipt.