

## Product Datasheet

### Rabbit IFN gamma protein (orb1215705)

<b>Description</b>	The Rabbit IFN gamma Biotinylated yeast-derived recombinant protein is not tagged and is naturally e...
<b>Target</b>	IFN gamma
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	QDTLTRETEHLKAYLKANTSDVANGGPLFLNILRNWKEESDNKIIQSQIVSFYFKL FDNLKDHEVIKKSMESEIKEDIFVKFFNSNLTKMDDFQNLTRISVDDRLVQRKAVS ELSNVLNFLSPKSNLKKRKRKRSQTLFRGRRASKY (144)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	16.9 kDa
<b>Entrez</b>	<a href="#">100008602</a>
<b>Expiration Date</b>	6 months from date of receipt.