

## Product Datasheet

### Canine Amphiregulin protein (orb1215712)

<b>Description</b>	The Canine Amphiregulin yeast-derived recombinant protein is not tagged and is naturally endotoxin-f...
<b>Target</b>	Amphiregulin
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	SVRVEQVVKP MKNKTESEKT SDKPKRKKKG GKSGKNRRNR KKKNPCDAEF QNFCIHGECK YIEHLEAVTC QCHQDYFGER CGEKSMK (87)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	10.1 kDa
<b>Entrez</b>	<a href="#">475176</a>
<b>Expiration Date</b>	6 months from date of receipt.