

## Product Datasheet

### Bovine IL-1 beta protein (orb1215715)

<b>Description</b>	The Bovine IL-1 beta Biotinylated yeast-derived recombinant protein is not tagged and is naturally e...
<b>Target</b>	IL-1 beta
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	APVQSIKCKLQDREQSLVLASPCVLKALHLLSQEMNREVVFCMSFVQGEERDN KIPVALGIKDKNLYLSCVKKGDPTLQLEEVDPKVYPKRNMEKRFVFKTEIKNTV EFESVLYPNWYISTSQIEERPVLGHFRGGQDITDFRMETLSP
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	17.7 kDa
<b>Entrez</b>	<a href="#">281251</a>
<b>Expiration Date</b>	6 months from date of receipt.