

Product Datasheet

Monkey IFN beta protein (orb1215726)

Description	The Squirrel Monkey IFN beta yeast-derived recombinant protein is not tagged and is naturally endoto...
Target	IFN beta
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	MSYNLLGFLQ RRSNFQSQKL LRQLNGRLEY CLKDRMDFEI PKEITQPQL QKEDVALIVY EMLQNIFTIF NRDFSRTGWN KTIVENFLAN VYQHIEHLKT ILKEKWEKKE FTREKSMSSL HLKRYYGRIL HYLMAKEYSH CAWTIVSAEI LRNFYFIDRL TGYL (164)
Purity	98%
Source	Yeast
MW	20.0 kDa
Entrez	101035681
Expiration Date	6 months from date of receipt.