

## Product Datasheet

### Mammal IFN Lambda protein (orb1215727)

<b>Description</b>	The Common Vampire Bat IFN lambda 3 yeast-derived recombinant protein is not tagged and is naturally...
<b>Target</b>	IFN Lambda
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	PPLLTSLMP SLVSPSINTD LPVGCTLVLM LMTTVLTRTA AVPVPTPLSA LPGARGCVVA QFKFLAPQDMKAFRRAKDTLEELLKPNRSCSSRPFRTRDLRQLQVWERPVAL EAELALTLKVLGSIANSTLGDILDQPLHMLRYIHTKLQACVPAQATAGPRPRGHLH HWLHRLQEASKKESTGCLEASVTFNLFRLTHDLQCVASGGICI (214)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	23.4 kDa
<b>Entrez</b>	<a href="#">112320026</a>
<b>Expiration Date</b>	6 months from date of receipt.