

## Product Datasheet

### Human IL-36 alpha protein (orb1215732)

<b>Description</b>	The Human IL-36 alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
<b>Target</b>	IL-36 alpha
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	KIDTPQQGSI QDINHRVWVL QDQTLIAVPR KDRMSPVTIA LISCRHVETL EKDRGNPIYL GLNGLNLCLM CAKVGDQPTL QLKEKDIMDL YNQPEPVKSF LFYHSQSGRN STFESVAFPG WFIAVSSEGG CPLILTQELG KANTTDFGLT MLF (153)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	17.1 kDa
<b>Entrez</b>	<b>27179</b>
<b>Expiration Date</b>	6 months from date of receipt.