

Product Datasheet

Bovine IL-17 protein (orb1215735)

Description	The Bovine IL-17A Biotinylated yeast-derived recombinant protein is not tagged and is naturally endo...
Target	IL-17
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	KAGVIIPQSP GCPPTEDKNF PQHVRVNLNI VNRSTNSRRP TDYHKRSTSP WTLHRNEDPE RYPSTVIWEAK CSHSGCINAE GKVDHMMNSV TIQQEILVLR RESQHCPHSF RLEKMLVAVG CTCVTPIVRH LA (132)
Purity	98%
Source	Yeast
MW	15.0 kDa
Entrez	282863
Expiration Date	6 months from date of receipt.