

Product Datasheet

Bovine IFN gamma protein (orb1215751)

Description The Bovine IFN gamma Biotinylated yeast-derived recombinant protein is not tagged and is naturally e...

Target IFN gamma

Form/Appearance Lyophilized

Storage -20°C

Note For research use only

Protein Sequence QGQFFREIEN LKEYFNASSP DVAKGGPLFS
EILKNWKDES DKKIIQSQIV SFYFKLFENL
KDNQVIQRSM DIIKQDMFQK FLNGSSEKLE
DFKKLIQIPV DDLQIQRKAI NELIKVMNDL
SPKSNLRKRK RSQNLFRGRR AST

Purity 98%

Source Yeast

MW 16.9 kDa

Entrez [281237](#)

Expiration Date 6 months from date of receipt.