

Product Datasheet

Porcine IFN alpha protein (orb1215752)

Description	The Swine IFN alpha 1 Biotinylated yeast-derived recombinant protein is not tagged and is naturally ...
Target	IFN alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	CDLPQTHSLAHTRALRLLAQMRRISPFSCLDHRRDFGSPHEAFGGNQVQKAQA MALVHEMLQQTFLFSTEGSAAAWNESLLHQFCTGLDQQLRDLEACVMQGAG LEGTPLLEEDSILAVRKYFHRLLTYLQEKSYSPCAWEIVRAEVMRSFSSSRNLQDR LRKKE (166)
Purity	98%
Source	Yeast
MW	19.0 kDa
Entrez	397686
Expiration Date	6 months from date of receipt.