

## Product Datasheet

### Equine IL-2 protein (orb1215755)

<b>Description</b>	The Equine IL-2 Biotinylated yeast-derived recombinant protein is not tagged and is naturally endoto...
<b>Target</b>	IL-2
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	APTSSSKRET QQQLKQLQMD LKLLLEGVNN NKNPKLSKML TFKINMPKKA TELKHLQCLE EELKPLEEML KNFLSKDIKE LMSNINVTVL GLKGSETRFT CEYDDEGTI VEFLNKWITF CQSIFSTMT
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	14.9 kDa
<b>Entrez</b>	<b>100034204</b>
<b>Expiration Date</b>	6 months from date of receipt.