

Product Datasheet

Ferret CCL16 protein (orb1215757)

Description	The Ferret CCL16 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	CCL16
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QSKSPASLNS SHTCCLKHHE KVMPRKMVVG YRRALNCYLP AIIFVTKKNR EICANPQNKW VQEYLRDPNI PLLPPRNLAQ ARSVRP (86)
Purity	98%
Source	Yeast
MW	9.9 kDa
Entrez	101674669
Expiration Date	6 months from date of receipt.