

Product Datasheet

Guinea pig IFN alpha 17 protein (orb1215759)

Description	The Guinea Pig IFN alpha 17 yeast-derived recombinant protein is not tagged and is naturally endotox...
Target	IFN alpha 17
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	CTLPQTHSQG NTRALILLEQ MRRISPFSCL KHRQDFGLAQ LDFEGKQVQK AQALSILHEM TRQTFNLFIS EDSSAAWNKS LLHSFCTGLH QQLNDLEACL TQETGLEEPP LMHEDSRLVV RKYFHRIAVY LKEKSYSPCA WEVVRAQIVR SLSASEKLAR KMQD (164)
Purity	98%
Source	Yeast
MW	18.9 kDa
Entrez	100717234
Expiration Date	6 months from date of receipt.