

Product Datasheet

Porcine CXCL10 protein (orb1215761)

Description	The Swine CXCL10 Biotinylated yeast-derived recombinant protein is not tagged and is naturally endot...
Target	CXCL10
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	PLSRTVRCTC IKISDRPVNP RSLEKLEMIP ASQSCPHVEI IATMKKNGEK RCLNPESKTI KNLLKAISKE RSKRSPRTQR EA
Purity	98%
Source	Yeast
MW	9.3 kDa
Entrez	494019
Expiration Date	6 months from date of receipt.