

Product Datasheet

Ferret CX3CL1 protein (orb1215766)

Description	The Ferret/Mink CX3CL1 yeast-derived recombinant protein is not tagged and is naturally endotoxin-fr...
Target	CX3CL1
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QHHGVKKCSV TCNKMTSEIP VALLDHYQRN QESCGKPAII LRTKKNRIFC ADPEAEWVRR AITHLDRQTS VPTRDG (76)
Purity	98%
Source	Yeast
MW	8.7 kDa
Entrez	101687656
Expiration Date	6 months from date of receipt.