

Product Datasheet

Monkey CTLA-4 protein (orb1215769)

Description	The Cynomolgus Monkey sCTLA-4 yeast-derived recombinant protein is not tagged and is naturally endot...
Target	CTLA-4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	KAMHVAQPAV VLANSRGIAS FVCEYASPGK ATEVRVTVLR QADSQVTEVC AATYMMGNEL TFLDDSICTG TSSGNQVNLT IQGLRAMDTG LYICKVELMY PPPYYMGIGN GTQIYVIDPE PCPDS (126)
Purity	98%
Source	Yeast
MW	13.5 kDa
Entrez	102115124
Expiration Date	6 months from date of receipt.