

Product Datasheet

Human B7-H2 protein (orb1215803)

Description	The Human B7-H2 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	B7-H2
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	DTQEKEVRAM VGSDVELSCA CPEGSRFDLN DVYVYWQTSE SKTVVTYHIP QNSSLENVDS RYRNRALMSP AGMLRGDFSL RLFNVTPQDE QKFHCLVLSQ SLGFQEVLSV EVTLHVAANF SVPVVSAPHS PSQDELTFTC TSINGYPRPN VYWINKTDNS LLDQALQNDT VFLNMRGLYD VVSVLRIART PSVNIGCCIE NVLLQQLTV GSQTGNDIGE RDKITENPVS TGEKNAAT (238)
Purity	98%
Source	Yeast
MW	26.4 kDa
Entrez	23308
Expiration Date	6 months from date of receipt.