

Product Datasheet

Human Frizzled-1 protein (orb1215809)

| | |
|-------------------------|---|
| Description | The Human Frizzled-1 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free... |
| Target | Frizzled-1 |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | QAAGQGPQGPGPGQQPPPPQQQSGQQYNGERGISVPDHGYCQPISIP LCT DIAYNQTIMP NLLGHTNQEDAGLEVHQFYPLVKVQCSAELKFFLC SMYAPVCTVL EQALPPCRSLCERARQGCEALMNKFGFQWPD TLKCEKFPVHGAGELCVGQNTS DKGTPTPSLLPEFWTSN (178) |
| Purity | 98% |
| Source | Yeast |
| MW | 19.4 kDa |
| Entrez | 8321 |
| Expiration Date | 6 months from date of receipt. |